Instructions for use

BIFURCATED NEEDLE

MD REF PB-02000-000-00

Roehling Medical Lancaster, LLC
44 Denver Road
Denver, PA 17517 USA
+1 (717) 335-3700

Medical Technology Promedt Consulting GmbH
Ernst-Heckel-Straße 7
66386 St. Ingbert, Germany.
ear@mt-procons.com

Indications: The **Roehling Medical Lancaster, LLC BIFURCATED ALLERGY SKIN TESTING NEEDLE** is a sterile, disposable, double pointed needle for in-vivo diagnostic skin testing using the puncture technique.
Description: The bifurcated needle is 6.6 cm long, with 0.5mm between the tines. Each needle is individually sterilized and sealed in 5-strip packs, available in packages of 100.

DIRECTIONS

- Wipe patient's skin with 70% ethanol or isopropyl alcohol and allow to air dry.
- Mark test sites with a pen, placing them 6cm apart for accurate wheal measurements.
- Place a drop of the test antigen next to each site.
- Hold the shaft of a sterile bifurcated needle between thumb and forefinger, and lightly puncture the skin at 90° angle. Penetrate the skin without drawing blood.
- Apply the positive and negative controls near the sites of the test allergens.

Because skin responses to mechanical and chemical stimulation can differ from patient to patient, we recommend using a positive control of 1mg/ml of histamine base or 5mg/ml of histamine phosphate, and a negative control of 50% glycerol saline, or whatever diluent is supplied with the allergenic extract.

- Grade the diagnostic site against the positive and negative sites. Record the readings 15 minutes after the puncture test and re-inspect the sites in 30 minutes. If the test sight has grown, discard the original reading. The positive control may determine whether a patient's reaction is impaired by certain drugs or for other reasons.

GRADING SYSTEM

Grade	Wheal (mm) after subtracting any wheal response to diluent control.
0	No Wheal
1+	< than half the histamine diameter
2+	half the size of histamine same size as histamine control
3+	same size as histamine control
4+	size of histamine control +2mm
5+	> than 2x histamine control

⚠️ PRECAUTIONS FOR USE

- Administer each allergy test with a separate bifurcated allergy skin testing needle and then discard the needle into a sharps container.
- Discard opened or damaged packaging.
- Report serious incidents to the manufacturer and relevant health authorities using the numbers listed above.
- If the needle is reused, the potential exists for the transmission of infectious agents which can result in serious infection or in extreme cases, death. Additionally with each reuse, the needle may become dull potentially resulting in inadequate skin penetration which may lead to incomplete delivery of expected dose. In addition, the packaging and/or instructions for use may not be present if the needle is reused which may result in improper use of the device leading to the possibility of the consequences previously stated.
- The bifurcated needle may be used with any commercially available allergenic extract for in-vivo diagnostic puncture skin testing to assess a patient's reactivity to a particular allergen. Roehling Medical Lancaster, LLC certifies that the bifurcated needle is sterile only when the seal of an individual unit is unbroken. Do not re-sterilize the needle, and take all necessary precautions associated with skin testing before using it. Federal law restricts this needle to sale by or on the order of a physician.

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ZWEIZINIGE NADEL

MD REF PB-02000-000-00

Roehling Medical Lancaster, LLC
44 Denver Road
Denver, PA 17517 USA
+1 (717) 335-3700

Medical Technology Promedt Consulting GmbH
Ernst-Heckel-Straße 7
66386 St. Ingbert, Germany.
ear@mt-procons.com

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